


# Bipolar Forceps

Exchangeable



These instrument can be disassembled into 4 pieces for easy cleaning

The electrode can be set in 4 positions ( every 90°)


Special setting of the handle, makes it possible to operate right or left handed

All elements of these Forceps, are 100 % compatible, they are very easy to repair ,at low cost

These instruments are completely manufactured in Germany

They have a very high quality performance

Because of our high quality, we give you a 1 year warranty.





- MM-25-2153 Forceps completed, with hook, working length: 33 cm
- MM-25-2154 Forceps completed, with hook, working length: 42 cm
- MM-25-2153-E Electrode only, with hook, working length: 33 cm
- MM-25-2154-E Electrode only, with hook, working length: 42 cm



- MM-25-2143 Forceps completed, without hook, working length: 33 cm
- MM-25-2144 Forceps completed, without hook, working length: 42 cm
- MM-25-2144-E Electrode only, without hook, working length: 33 cm
- MM-25-2144-E Electrode only, without hook, working length: 42 cm



- MM-25-2163 Forceps completed, small jaw, working length: 33 cm
- MM-25-2164 Forceps completed, small jaw, working length: 42 cm
- MM-25-2163-E Electrode only, small jaw, working length: 33 cm
- MM-25-2164-E Electrode only, small jaw, working length: 42 cm



- MM-25-2173 Forceps completed, dolphin jaw, working length: 33 cm
- MM-25-2174 Forceps completed, dolphin jaw, working length: 42 cm
- MM-25-2173-E Electrode only, dolphin jaw, working length: 33 cm
- MM-25-2174-E Electrode only, dolphin jaw, working length: 42 cm



- MM-25-2183 Forceps completed, wide jaw, working length: 33 cm
- MM-25-2184 Forceps completed, wide jaw, working length: 42 mm
- MM-25-2183-E Electrode only, wide jaw, working length: 33 cm
- MM-25-2184-E Electrode only, wide jaw, working length: 42 cm

## Bipolar Cables



MM-25-2100 Connexion to Martin and Berchtold units, length: 3mm

MM-25-2101 Connexion to Martin and Berchtold units, length: 5mm



MM-25-2102 Connexion to Erbe and Select units, length: 3mm

MM-25-2103 Connexion to Erbe and Select units, length: 5mm



MM-25-2104 Connexion to US units, length: 3mm

MM-25-2105 Connexion to US units, length: 5mm

## Spare Parts

MM-25-2100-8 Rubber tip

MM-25-2103-70 Sheath, plastic, 33 cm

MM-25-2104-70 Sheath, plastic, 42 cm

MM-25-2105-70 Sheath, metal, 33 cm

MM-25-2106-70 Sheath, metal, 42 cm