


Fluid-Light-Cable



Our Fluid-Light-Cable gives you 100 % more light and a higher color temperature.

No transmission errors with broken fibers.

No fiber glass bundle, but a special fluid medium in a high reflecting and a flexible tube.

Fluid-Light-Cables are only soakable. Do not sterilize them with steam or gas.

Fluid-Light-Cables

MM-25-Z-5181 Fluid-Light-Cable, 5mm diam., Lenght: 180 cm
MM-25-Z-5281 Fluid-Light-Cable, 5mm diam., Lenght: 230 cm

Adaptor



MM-25-Z-5100 Adaptor Projector Storz/Wisap
MM-25-Z-5101 Adaptor Projector Wolf/Dyonics
MM-25-Z-5102 Adaptor Projector ACMI
MM-25-Z-5103 Adaptor Projector Olympus



MM-25-Z-5110 Adaptor Instrument Storz/Olympus/Wisap
MM-25-Z-5111 Adaptor Instrument Wolf/Dyonics
MM-25-Z-5110 Adaptor Instrument ACMI