



# Instruments for Laparoscopy

Rotatable and exchangeable





## EXAMPLE




### Features "MULTIFUNCTION"

- High Precision & Advanced Technology
  - Weight: 100 - 120g
  - Full Autoclavable
  - Interchangeable jaws
    - 3 parts
- Easy to assemble and disassemble
- Shaft rotates in 15° increments
  - Insulated
- Lower costs are obtained by using multiple jaws with one handle
- Absolutely Compatibility between 5 and 10mm Instruments







**P9-070**  
Handle without ratchet,  
insulated for 5mm and 10mm




**P9-071**  
Handle with ratchet,  
insulated for 5mm and 10mm




**P9-072**  
Handle with spring ratchet,  
insulated for 5mm and 10mm



**P9-073**  
Spare part  
only for 5mm sheath (P9-074)




**P9-074**  
5mm shaft for dismantable instruments









**P9-075**  
10mm shaft for dismantable instruments










**MM-P9-076**  
Rotating cap for Handle  
MM-P9-070 MM-P9-071  
MM-P9-072









Length 330mm and 450mm







| Grasping & Dissecting Forceps 5mm  | Part - No.:  |
|--|--|
|  <p>Serrated, cupped jaws, 5mm</p>                            | <p>MM-P9-125-33<br/>MM-P9-125-45</p> <p>Double-ation</p> |
|  <p>Crps serrated artrumatic, 5mm</p>                         | <p>MM-P9-126-33<br/>MM-P9-126-45</p> <p>Double-ation</p> |
|  <p>Dorsey, fenestrated, 5mm</p>                             | <p>MM-P9-127-33<br/>MM-P9-127-45</p> <p>Double-ation</p> |
|  <p>Allis, 5mm</p>  | <p>MM-P9-128-33<br/>MM-P9-128-45</p> <p>Double-ation</p> |
|  <p>Dorsey Delicate, serrated and fenestrated jaws, 5mm</p> | <p>MM-P9-129-33<br/>MM-P9-129-45</p> <p>Double-ation</p> |
|  <p>Atraumatic, serrated, 5mm</p>                           | <p>MM-P9-130-33<br/>MM-P9-130-45</p> <p>Single-ation</p> |


| Grasping & Dissecting Forceps 5mm   | Part - No.:  |
|---|--|
|  <p>Serrated, conical jaws with cup, 5mm</p>            | <p>MM-P9-131-33<br/>MM-P9-131-45</p> <p>Double-ation</p> |
|  <p>Mc Keman, cross and right angles serrated, 5mm</p>  | <p>MM-P9-132-33<br/>MM-P9-132-45</p> <p>Double-ation</p> |
|  <p>Reddick-Olsen serrated, 5mm</p>                    | <p>MM-P9-133-33<br/>MM-P9-133-45</p> <p>Double-ation</p> |
|  <p>Alligator forceps, serrated jaws 5mm</p>          | <p>MM-P9-141-33<br/>MM-P9-141-45</p> <p>Double-ation</p> |
|  <p>Single-ation</p>                                  | <p>MM-P9-142-33<br/>MM-P9-142-45</p> <p>Single-ation</p> |
|  <p>Cobra, serrated jaws with delicate teeth, 5mm</p> | <p>MM-P9-143-33<br/>MM-P9-143-45</p> <p>Single-ation</p> |









Length 330mm and 450mm







| Grasping & Dissecting Forceps 5mm   | Part - No.:  |
|---|--|
|  <p>Manhes, delicate serrated, 5mm</p>   | <p>MM-P9-144-33<br/>MM-P9-144-45</p> <p>Single-action</p>  |
|  <p>Smooth, 5mm</p>                      | <p>MM-P9-145-33<br/>MM-P9-145-45</p> <p>Double -action</p> |
|  <p>Serrated with teeth, 5mm</p>        | <p>MM-P9-146-33<br/>MM-P9-146-45</p> <p>Double -action</p> |
|  <p>Manhes 1x2 teeth serrated, 5mm</p> | <p>MM-P9-147-33<br/>MM-P9-147-45</p> <p>Single-action</p>  |
|  <p>Babcock double -action, 5mm</p>    | <p>MM-P9-148-33<br/>MM-P9-148-45</p> <p>Double -action</p> |
|  <p>Alligator forceps, 5mm</p>         | <p>MM-P9-149-33<br/>MM-P9-149-45</p> <p>Double -action</p> |


| Grasping & Dissecting Forceps 5mm   | Part - No.:  |
|---|--|
|  <p>Wave grasper, 5mm</p>                       | <p>MM-P9-150-33<br/>MM-P9-150-45</p> <p>Single-action</p>  |
|  <p>Delicate serrated long forceps 5mm</p>      | <p>MM-P9-151-33<br/>MM-P9-151-45</p> <p>Double -action</p> |
|  <p>Alligator, strong teeth forceps 5mm</p>    | <p>MM-P9-152-33<br/>MM-P9-152-45</p> <p>Single-action</p>  |
|  <p>Right angel, 5mm</p>                      | <p>MM-P9-153-33<br/>MM-P9-153-45</p> <p>Double -action</p> |
|  <p>Fenestrated, 5mm</p>                      | <p>MM-P9-154-33<br/>MM-P9-154-45</p> <p>Single-action</p>  |
|  <p>Maryland, serrated, curved left, 10mm</p> | <p>MM-P9-155-33<br/>MM-P9-155-45</p> <p>Double -action</p> |









Length 330mm and 450mm



| Grasping & Dissecting Forceps 5mm / 10mm  | Part - No.: |
|---|-------------|
|  <p>MM-P9-156-33<br/>MM-P9-156-45</p> <p>Mixer 90° angled, cross-serrated, 10mm Double-action</p>      |             |
|  <p>MM-P9-157-33<br/>MM-P9-157-45</p> <p>Maryland, serrated curved left, 5mm Double-action</p>         |             |
|  <p>MM-P9-158-33<br/>MM-P9-158-45</p> <p>De Bakey middle 40mm atraum, serrated 10mm Double-action</p> |             |
|  <p>MM-P9-159-33<br/>MM-P9-159-45</p> <p>Glassman, atraumatic serrated 10mm Double-action</p>        |             |
|  <p>MM-P9-160-33<br/>MM-P9-160-45</p> <p>Mixer, left curved, 10mm Double-action</p>                  |             |
|  <p>MM-P9-161-33<br/>MM-P9-161-45</p> <p>Overholt (Kelly) left curved, 10mm Double-action</p>        |             |


| Grasping & Dissecting Forceps 5mm / 10mm   | Part - No.: |
|--|-------------|
|  <p>MM-P9-162-33<br/>MM-P9-162-45</p> <p>De Bakey middle 30mm atraum, serrated 10mm Double-action</p>      |             |
|  <p>MM-P9-163-33<br/>MM-P9-163-45</p> <p>De Bakey middle 50mm atraum, serrated 10mm Double-action</p>      |             |
|  <p>MM-P9-164-33<br/>MM-P9-164-45</p> <p>Babcock, 10mm Double-action</p>                                  |             |
|  <p>MM-P9-165-33<br/>MM-P9-165-45</p> <p>Wave Grasper, 10mm Double-action</p>                            |             |
|  <p>MM-P9-167-33<br/>MM-P9-167-45</p> <p>Claw Grasper, large jaws, with 2x3 teeth 10mm Single-action</p> |             |
|  <p>MM-P9-168-33<br/>MM-P9-168-45</p> <p>Dorsey, fenestrated, 10mm Double-action</p>                     |             |





Length 330mm and 450mm



| Scissors 5mm / 10mm   | Part - No.: |
|---|-------------|
|  <p>MM-P9-169-33<br/>MM-P9-169-45</p> <p>Micro Scissors,<br/>jaws straght, 5mm</p> <p><b>Single-action</b></p>                                 |             |
|  <p>MM-P9-170-33<br/>MM-P9-170-45</p> <p>Hook Scissors,<br/>5mm</p> <p><b>Single-action</b></p>  |             |
|  <p>MM-P9-171-33<br/>MM-P9-171-45</p> <p>Metzenb.Scissors,<br/>curved left, 5mm</p> <p><b>Double-action</b></p>                               |             |
|  <p>MM-P9-172-33<br/>MM-P9-172-45</p> <p>Scissors Mayo Model one<br/>blade serratted,<br/>straght jaws, 10mm</p> <p><b>Single-action</b></p> |             |
|  <p>MM-P9-173-33<br/>MM-P9-173-45</p> <p>Micro Scissors,<br/>jaws curved left,<br/>5mm</p> <p><b>Single-action</b></p>                       |             |
|  <p>MM-P9-174-33<br/>MM-P9-174-45</p> <p>Scissors, short, serrated<br/>jaws, 5mm</p> <p><b>Single-action</b></p>                             |             |

| Scissors 5mm / 10mm  | Part - No.: |
|--|-------------|
|  <p>MM-P9-175-33<br/>MM-P9-175-45</p> <p>metzenb.Scissors,<br/>straght jaws, 5mm</p> <p><b>Single-action</b></p>                                 |             |
|  <p>MM-P9-176-33<br/>MM-P9-176-45</p> <p>Scissors Mayo Model one<br/>blade serratted,<br/>left curved, jaws 10mm</p> <p><b>Single-action</b></p> |             |
|  |             |
|  |             |
|  |             |
|  |             |




Lenght 330mm and 450mm




| Biopsy Forceps<br>5mm / 10mm  | Part - No.:   |
|---|---|
|  | <p>MM-P9-177-33<br/>MM-P9-177-45</p> <p>Serrated cupped jaws,<br/>pin, 5mm</p> <p>Double-action</p>     |
|  | <p>MM-P9-178-33<br/>MM-P9-178-45</p> <p>Sharp cupped jaws without<br/>pin, 5mm</p> <p>Double-action</p> |
|   |   |
|   |   |
|   |   |
|   |   |



| Biopsy Forceps<br>5mm / 10mm   | Part - No.:  |
|--|--|
|  | <p>MM-P9-179-33<br/>MM-P9-179-45</p> <p>Sharp cupped jaws with<br/>pin, 5mm</p> <p>Double-action</p>     |
|  | <p>MM-P9-180-33<br/>MM-P9-180-45</p> <p>Sharp cupped jaws without<br/>pin, 10mm</p> <p>Single-action</p> |
|  |  |
|  |  |
|  |  |
|  |  |









Length 330mm and 450mm


| Biopsy Forceps<br>5mm / 10mm:  | Part - N.:                   |
|--|------------------------------|
|   | MM-P9-181-33<br>MM-P9-181-45 |
| Biopsy Forceps cutting,<br>5mm   | Single action                |
|   | MM-P9-182-33<br>MM-P9-182-45 |
| Hook Biopsy Forceps,<br>5mm  | Single action                |
|  | MM-P9-183-33<br>MM-P9-183-45 |
| Biopsy Forceps cutting,<br>5mm   | Single action                |
|  |                              |
|  |                              |
|  |                              |

| Biopsy Forceps<br>5mm / 10mm:  | Part - N.:                   |
|--|------------------------------|
|  | MM-P9-185-33<br>MM-P9-185-45 |
| Biopsy Forceps cutting with<br>pin, 5mm  | Single action                |
|  | MM-P9-186-33<br>MM-P9-186-45 |
| Biopsy Forceps cutting,<br>5mm   | Single action                |
|  |                              |
|  |                              |
|  |                              |

**Special Needleholder with ergonomic handle**
















| Part - No.:   |  |
|---|--|
|  <p>MM-P9-236</p> <p>Needleholder, straight, with tungsten carbide jaws, 5mm</p>       |  |
|  <p>MM-P9-237</p> <p>Needleholder, curved left, with tungsten carbide jaws, 5mm</p>  |  |
|  <p>MM-P9-238</p> <p>Needleholder, curved right, with tungsten carbide jaws, 5mm</p> |  |



|   |  |   |  |  |
|---|--|---|--|--|
| <p>Traumatic grasping Forceps, single action jaw<br/>φ 5 mm</p>  <p>MM-21-0100-33/45</p>       | <p>Traumatic grasping Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-010-33/45</p>           | <p>Atraumatic grasping Forceps, single action jaw<br/>φ 5 mm</p>  <p>MM-21-0102-33/45</p>        | <p>Atraumatic grasping Forceps, single action jaw<br/>φ 5 mm</p>  <p>MM-21-0103-33/45</p>        | <p>Atraumatic grasping Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0104-33/45</p>         |
| <p>Kelly grasping Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0105-33/45</p>           | <p>Grasping Forceps, conical, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0106-33/45</p>           | <p>Atraumatic grasping Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0107-33/45</p>        | <p>Grasping Forceps, fine serrated double action jaw<br/>φ 5 mm</p>  <p>MM-21-0108-33/45</p>     | <p>Grasping Forceps, heavy serrated single action jaws<br/>φ 3 mm</p>  <p>MM-21-0109-33/45</p>    |
| <p>Grasping Forceps, heavy serrated double action jaw<br/>φ 5 mm</p>  <p>MM-21-0110-33/45</p>  | <p>Atraumatic grasping Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0111-33/45</p>         | <p>Biopsy Forceps, single action jaw<br/>φ 3 mm</p>  <p>MM-21-0112-33/45</p>                     | <p>Biopsy Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0113-33/45</p>                     | <p>Biopsy Forceps, single action jaw<br/>φ 3 mm</p>  <p>MM-21-0114-33/45</p>                      |
| <p>Biopsy Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0115-33/45</p>                   | <p>Large spoon Forceps, single action jaw<br/>φ 10 mm</p>  <p>MM-21-0116-33/45</p>                | <p>Biopsy Forceps, long jaw single action jaw<br/>φ 5 mm</p>  <p>MM-21-0117-33/45</p>            | <p>Biopsy Forceps, with thorn long jaw single action jaw<br/>φ 5 mm</p>  <p>MM-21-0118-33/45</p> | <p>Biopsy Forceps, with teeth single action jaw<br/>φ 3 mm</p>  <p>MM-21-0119-33/45</p>           |
| <p>Biopsy Forceps, with teeth double action jaw<br/>φ 5 mm</p>  <p>MM-21-0120-33/45</p>        | <p>Fenestrated Forceps, heavy serrated, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0121-33/45</p> | <p>Fenestrated Forceps, fine serrated, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0122-33/45</p> | <p>Grasping Forceps, serrated spoon single action jaw<br/>φ 5 mm</p>  <p>MM-21-0123-33/45</p>    | <p>Grasping Forceps, serrated 2x3 teeth single action jaw<br/>φ 5 mm</p>  <p>MM-21-0124-33/45</p> |
| <p>Maryland Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0125-33/45</p>                | <p>Mixer Forceps, double action jaw<br/>φ 5 mm</p>  <p>MM-21-0126-33/45</p>                      | <p>Hook Scissor single action jaw<br/>φ 5 mm</p>  <p>MM-21-0127-33/45</p>                       | <p>Scissors single action jaw<br/>φ 5 mm</p>  <p>MM-21-0128-33/45</p>                           | <p>Metzenbaum Scissors Straight double action jaw<br/>φ 5 mm</p>  <p>MM-21-0129-33/45</p>        |
| <p>Metzenbaum Scissors left curved double action jaw<br/>φ 5 mm</p>  <p>MM-21-0130-33/45</p> | <p>Micro Scissors single action jaw<br/>φ 5 mm</p>  <p>MM-21-0131-33/45</p>                     | <p>Large Operating Scissors φ 10 mm</p>  <p>MM-21-0132-33/45</p>                               | <p>Babcock Clamp double action jaw<br/>φ 5 mm</p>  <p>MM-21-0133-33/45</p>                    |  |
| <p>Babcock spezial jaw double action jaw<br/>φ 10 mm</p>  <p>MM-21-0134-33/45</p>            | <p>Duval Clamp double action jaw<br/>φ 10 mm</p>  <p>MM-21-0135-33/45</p>                       |   | <p>Allis Clamp double action jaw<br/>φ 5 mm</p>  <p>MM-21-0136-33/45</p>                       |  |
| <p>DeBaKey Clamp double action jaw<br/>φ 5 mm</p>  <p>MM-21-0137-33/45</p>                   |  | <p>Pennington Clamp double action jaw<br/>φ 5 mm</p>  <p>MM-21-0138-33/45</p>                 |  | <p>Biopsy Punch single action jaw<br/>φ 5 mm</p>  <p>MM-21-0139-33/45</p>                       |
| <p>Biopsy Punch with thorn single action jaw<br/>φ 5 mm</p>  <p>MM-21-0140-33/45</p>         | <p>Hook Punch single action jaw<br/>φ 5 mm</p>  <p>MM-21-0141-33/45</p>                         | <p>Grasping Forceps serrated single action jaw<br/>φ 3 mm</p>  <p>MM-21-0142-33/45</p>         | <p>Grasping Forceps serrated double action jaw<br/>φ 5 mm</p>  <p>MM-21-0143-33/45</p>         | <p>Tapered grasping Forceps single action jaw<br/>φ 3 mm</p>  <p>MM-21-0144-33/45</p>           |

|  |   |  |   |  |
|--|---|--|---|--|
| <p>Tapered Forceps<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0145-33/45</p>                            | <p>Delphin nose Forceps<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0146-33/45</p>                          | <p>Fundus Forceps<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0147-33/45</p>                         | <p>Claw Forceps<br/>single action jaw<br/>φ 5 mm</p>  <p>MM-21-0148-33/45</p>                                | <p>Atraumatic Clamp<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0149-33/45</p>                                      |
| <p>Atraumatic Clamp<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-0150-33/45</p>                          |   | <p>Dorsey Clamp<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0151-33/45</p>                           | <p>Needle holde<br/>φ 5 mm</p>  <p>MM-21-0152-33/45</p>   | <p>Micro Forceps<br/>φ 5 mm</p>  <p>MM-21-0153-33/45</p>  |
| <p>Dessector, serrated<br/>φ 5 mm</p>  <p>MM-21-0154-33/45</p>  | <p>Dessector, serrated<br/>strong curve<br/>φ 5 mm</p>  <p>MM-21-0155-33/45</p>                                | <p>Grasping Forceps, blunt<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0156-33/45</p>                | <p>Mixer grasping Forceps,<br/>45° angle double<br/>action jaw φ 10 mm</p>  <p>MM-21-0157-33/45</p>         | <p>Fenestrated Atraumatic<br/>Clamp double action<br/>jaw φ 5 mm</p>  <p>MM-21-0158-33/45</p>                         |
| <p>Atraumatic Grasping Forceps,<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0159-33/45</p>               | <p>Delphin nose, short<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0160-33/45</p>                           | <p>Delphin nose, long<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0161-33/45</p>                     | <p>Grasping Forceps, Alligator<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0162-33/45</p>               | <p>Metzenbaum Scissors<br/>straight<br/>φ 5 mm</p>  <p>MM-21-0163-33/45</p>   |
| <p>Metzenbaum Scissors<br/>curved<br/>φ 5 mm</p>  <p>MM-21-0164-33/45</p>                                   | <p>Metzenbaum Scissors<br/>straight serrated<br/>φ 5 mm</p>  <p>MM-21-0165-33/45</p>                           | <p>Metzenbaum Scissors<br/>curved serrated<br/>φ 5 mm</p>  <p>MM-21-0166-33/45</p>                      | <p>DeBakey Forceps<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0167-33/45</p>                            | <p>Kocher Forceps 40 mm jaw<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-0168-33/45</p>                            |
| <p>Carbide needle holder<br/>10 mm jaw<br/>φ 5 mm</p>  <p>MM-21-0169-33/45</p>                             | <p>Carbide needle holder<br/>8 mm jaw<br/>φ 5 mm</p>  <p>MM-21-0170-33/45</p>                                 | <p>Grasping Forceps,<br/>single action jaw<br/>φ 10 mm</p>  <p>MM-21-0171-33/45</p>                    | <p>Claw Forceps,<br/>single action jaw<br/>φ 10 mm</p>  <p>MM-21-0172-33/45</p>                            | <p>Operating Scissors<br/>φ 10 mm</p>  <p>MM-21-0173-33/45</p>  |
| <p>Metzenbaum Scissors<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-0174-33/45</p>                     | <p>Peritoneal Scissors<br/>single action jaw<br/>φ 5 mm</p>  <p>MM-21-0175-33/45</p>                         | <p>Metzenbaum Scissors<br/>serrated<br/>φ 5 mm</p>  <p>MM-21-0176-33/45</p>                           | <p>Atraumatic Grasping Forceps,<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0177-33/45</p>             | <p>Crocodile Grasping Forceps,<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0178-33/45</p>                        |
| <p>Hook electrode 260 mm<br/>φ 5 mm</p>  <p>MM-21-0179-33/45</p>  | <p>J-Hook electrode 260 mm<br/>φ 5 mm</p>  <p>MM-21-0180-33/45</p>   | <p>Allis Atraumatic<br/>Grasping Forceps<br/>φ 5 mm</p>  <p>MM-21-0181-33/45</p>                      | <p>Hook electrode<br/>with angled shaft<br/>φ 5 mm</p>  <p>MM-21-0182-33/45</p>                          | <p>Electrode<br/>with angled shaft<br/>φ 5 mm</p>  <p>MM-21-0183-33/45</p>  |
| <p>J-Hook electrode<br/>with angled shaft<br/>φ 5 mm</p>  <p>MM-21-0184-33/45</p>                         | <p>Metzenbaum Scissors<br/>angled shaft single action jaw<br/>φ 5 mm</p>  <p>MM-21-0185-33/45</p>            | <p>Peritoneal Scissors<br/>angled shaft single action jaw<br/>φ 5 mm</p>  <p>MM-21-0186-33/45</p>     | <p>Micro preparation Scissors<br/>angled shaft single action jaw<br/>φ 5 mm</p>  <p>MM-21-0187-33/45</p>  | <p>Cobra Grasping Forceps<br/>with teeth angled shaft<br/>double action jaw<br/>φ 5 mm</p>  <p>MM-21-0188-33/45</p> |
| <p>Atraumatic Grasping Forceps<br/>angled shaft double action jaw<br/>φ 5 mm</p>  <p>MM-21-0189-33/45</p> | <p>Babcock/allis grasping Forceps<br/>angled shaft double action jaw<br/>φ 5 mm</p>  <p>MM-21-0190-33/45</p> | <p>Duval Grasping Forceps<br/>angled shaft double action jaw<br/>φ 10 mm</p>  <p>MM-21-0191-33/45</p> | <p>Alligator Grasping Forceps<br/>angled shaft double action jaw<br/>φ 5 mm</p>  <p>MM-21-0192-33/45</p> | <p>Atraumatic Grasping Forceps<br/>angled shaft double action jaw<br/>φ 5 mm</p>  <p>MM-21-0193-33/45</p>           |
| <p>Biopsy Punch<br/>single action jaw<br/>φ 3 mm</p>  <p>MM-21-0194-33/45</p>                             | <p>Kocher Forceps 60 mm jaw<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-0195-33/45</p>                   | <p>Hook Scissors<br/>single action jaw<br/>φ 3mm</p>  <p>MM-21-0196-33/45</p>                         | <p>Metzenbaum Scissors left curved<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-0197-33/45</p>         | <p>Atraumatic grasping Forceps,<br/>double action jaw<br/>φ 3 mm</p>  <p>MM-21-0198-33/45</p>                       |

|  |  |   |   |  |
|--|--|---|---|--|
| <p>Alligator Forceps,<br/>double action jaw<br/>φ 3 mm</p>  <p>MM-21-0199-33/45</p> |  <p>MM-21-0200-33/45</p>  | <p>Grasping Forceps,serrated,<br/>2x3 teeth single action jaw<br/>φ 10 mm</p>  <p>MM-21-0201-33/45</p> | <p>Maryland Forceps,<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-0202-33/45</p> | <p>Mixer Forceps,<br/>double action jaw<br/>φ 10 mm</p>  <p>MM-21-01203-33/45</p> |
| <p>Babcock Clamp double action jaw<br/>φ 10 mm</p>  <p>MM-21-0204-33/45</p>         |  | <p>Allis Clamp double action jaw<br/>φ 10 mm</p>  <p>MM-2-0205-33/45</p>                              |   |  |
| <p>DeBakey Clamp double action jaw<br/>φ 10 mm</p>  <p>MM-2-206-33/45</p>           | <p>Pennington Clamp double action jaw<br/>φ 10 mm</p>  <p>MM-21-0207-33/45</p> |   | <p>Needle holder<br/>φ 3 mm</p>  <p>MM-21-0208-33/45</p>                           |  |
| <p>Claw Forceps,single action jaw<br/>φ 10 mm</p>  <p>MM-21-0209-33/45</p>          | <p>Biopsy Punch<br/>single action jaw<br/>φ 3 mm</p>  <p>MM-21-0210-33/45</p>   | <p>Fenestrated Forceps,fine<br/>serrated,double action jaw</p>  <p>MM-21-0211-37/45</p>               |   |  |
|  |  |   |   |  |
|  |  |   |   |  |
|  |  |   |   |  |
|  |  |   |   |  |
|  |  |   |   |  |
|  |  |   |   |  |
|  |  |   |   |  |