


# Urology-Gynecology



## Resectoscopes



MM-25-U-2700-P Working element, passive cutting by spring action, closed handle electrode is inside the sheath in rest position.

MM-25-U-2700-A Working element, active cutting by finger action, closed handle electrode is outside the sheath in rest position.



MM-25-U-2701-P Working element, passive cutting by spring action, open handle electrode is inside the sheath in rest position.

MM-25-U-2701-A Working element, active cutting by finger action, open handle electrode is outside the sheath in rest position.



MM-25-U-2700 HF-Connection for working elements, length: 300 cm



Resectoscope Sheath for continuous flow, consisting of 3 pieces:

- MM-25-U-2702-A Outer sheath, 27ch. with perforation
- MM-25-U-2702-I Inner sheath, 24 ch. with ceramic isolation
- MM-25-U-2702-O Standard abturator



- MM-25-U-2702-B Deflecting abturator, for 24 ch. sheaths (MM-25-U-2702)
- MM-25-U-2703-B Deflecting abturator, for 27 ch. sheaths (MM-25-U-2703)

- MM-25-U-2702-S Visual abturator, for use with 24 ch. sheaths
- MM-25-U-2703-S Visual abturator, for use with 27 ch. sheaths



- MM-25-U-2702-Z Res., cope sheath with central valve, 24 ch. w.stand. abturator (MM-25-U-2702-O)
- MM-25-U-2703-Z Res., cope sheath with central valve, 27 ch. w.stand. abturator (MM-25-U-2703-O)

## Electrodes



MM-25-U-2301 Vaporisation electrode, 24 ch



MM-25-U-2302 Roller.ball electrode, 24 ch



MM-25-U-2303 Knife electrode, 24 ch



MM-25-U-2304 Cutting-loop, straight,24 ch



MM-25-U-2305 Cutting-loop,angled, 24 ch



MM-25-U-2306 Conical electrode,24 ch

## Cystoscope-Urethroscopes




- MM-25-U-2619 Urethroscope sheath, 17 ch., with abturator
- MM-25-U-2621 Urethroscope sheath, 19 ch., with abturator
- MM-25-U-2623 Urethroscope sheath 21 ch., with abturator
- Urethroscope sheath, 23 ch., with abturator



MM-25-U-0100  
Telescope bridge



MM-25-U-0101  
Telescope bridge with  
1 Instrument channe




MM-25-U-0102  
Telescope bridge with  
2 Instrument channels



- MM-25-U-2401 Catheter deflecting Mechanism with 1 Instrument channel with ratchet
- MM-25-U-2402 Catheter deflecting Mechanism with 2 Instrument channel with ratchet

## Optical Instruments



MM-25-U-2501 Biopsy Forceps, double action jaws



MM-25-U-2502 Grasping Forceps, double action jaws



MM-25-U-2503 Biopsy punch, through-cutting double action jaws



MM-25-U-2504 Hook scissors, double action jaws



MM-25-U-2505 Biopsy Forceps, down, double action jaws



MM-25-U-2506 Stone Punch for use with 25 ch. sheaths



MM-25-U-2850 Operating sheath 7 mm  $\Phi$  with 2 stop cocks rotatable, with distance holder and 1 instrument channel


MM-25-U-2840 Diagnostic sheath 5 mm  $\Phi$  graduated 1 stop cocks fixed, with distance holder and obturator



MM-25-U-2810 Sheath, with continuous flow, 5 mm  
MM-25-U-2820 Sheath, with 1 Instrument channel, 7 ch. 7 mm  
MM-25-U-2830 Sheath, with 2 Instrument channel, 7 ch. 8 mm

# Flexible Instruments

7 ch., Working length: 40 cm



MM-25-U-29230 Biopsy forceps, double action jaws



MM-25-U-29231 Grasping forceps, double action jaws



MM-25-U-29232 Scissors, single action jaws

# Evacuator



MM-25-U-29330 Evacuator complete

MM-25-U-29331 Glass for MM-25-U-29330

MM-25-U-29332 Rubber ball for MM 25-U-29330

MM-25-U-29333 Rubber hose for MM-25-U-29330

MM-25-U-29334 Connection for MM-25-U-29330